

Enhancement of β -Alanine Biosynthesis in *Escherichia coli* Based on Multivariate Modular Metabolic Engineering

Jian Xu, Li Zhou, Zhemin Zhou

Abstract

β -alanine is widely used as an intermediate in industrial production. However, the low production of microbial cell factories limits its further application. Here, to improve the biosynthesis production of β -alanine in *Escherichia coli*, multivariate modular metabolic engineering was recruited to manipulate the β -alanine biosynthesis pathway through keeping the balance of metabolic flux among the whole metabolic network. The β -alanine biosynthesis pathway was separated into three modules: the β -alanine biosynthesis module, TCA module, and glycolysis module. Global regulation was performed throughout the entire β -alanine biosynthesis pathway rationally and systematically by optimizing metabolic flux, overcoming metabolic bottlenecks and weakening branch pathways. As a result, metabolic flux was channeled in the direction of β -alanine biosynthesis without huge metabolic burden, and 37.9 g/L β -alanine was generated by engineered *Escherichia coli* strain B0016-07 in fed-batch fermentation. This study was meaningful to the synthetic biology of β -alanine industrial production.



Winpact Fermentation System

Introduction

E. coli was rationally and systematically engineered for the production of β -alanine. The MME strategy was applied to channel the metabolic flux in the direction of β -alanine biosynthesis with the aim of maintaining the balance of the intracellular metabolic network. The β -alanine biosynthesis pathway was separated into three modules: the β -alanine biosynthesis module, TCA module, and glycolysis module. Global regulation was performed throughout the entire β -alanine biosynthesis pathway rationally and systematically; as a result, 37.9 g/L β -alanine was generated in fed-batch fermentation.

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Room 2519, No. 306 Taopu Road,
Putuo District, Shanghai, China
Postal code: 200333
Cell: +86-13621808204
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D.No.7-143, 2nd Floor,
St.No.2, Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
T/ +91-40-27001586

Materials and Methods

Fed-batch fermentation was carried out at 37 °C in a 5 L Winpack Major Science bioreactor (working volume 2 L). DO was maintained above 45% by adjusting airflow (2–10 L/min) and agitation (200–900 rpm). Feeding was automatically regulated based on growth rate (μ) and glycerol consumption rate (qGly).

M9Y medium per liter contained: 1 g NH_4Cl , 0.5 g NaCl , 2 g yeast extract, 3 g KH_2PO_4 , 5 g glycerol, 5 mL metal solution, 6 g Na_2HPO_4 , and 13.21 g $(\text{NH}_4)_2\text{SO}_4$.

Feed medium per liter contained: 3.67 g MgSO_4 , 4 g yeast extract, 4 g tryptone, 100 g $(\text{NH}_4)_2\text{SO}_4$, and 650 g glycerol.

Results

IFed-batch fermentation was performed to provide a highly controlled environment, allowing the engineered *E. coli* strain to grow efficiently and accumulate a high concentration of the target product, β -alanine.

References

Evaluation of Metabolic Engineering Strategies on 2-Ketoisovalerate Production by *Escherichia coli*

<https://doi.org/10.3390/biology101010>

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Room 2519, No. 306 Taopu Road,
Putuo District, Shanghai, China
Postal code: 200333
Cell: +86-13621808204
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D.No.7-143, 2nd Floor,
St.No.2,Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
T/ +91-40-27001586

major
science
www.majorsci.com
info@majorsci.com