



2018 BROCHURE



Power Supply



New Products



BluView Transilluminator

MBE-300

MBE-200A



Major Science BluView Transilluminator uses the harmless blue LED lights to replace the aggressive UV lights, and allows you to directly view the experiment result without wearing any UV protection equipment.



* For more information, please refer to p20

Mini Horizontal Gel

MJ-105A



The MJ-105A horizontal electrophoresis system offers many advantages for nucleic acid separation. Agarose gels are convenient to cast and samples are easy to load with this recognizable mini horizontal cell.

* For more information, please refer to p13

MS Mini Incubator



MO-MINI



Major Science's Mini Incubator is designed for personal use and small laboratories, saving much of space. The unit features a broad temperature range to meet a variety of microbiology or hematology applications.

* For more information, please refer to p24

Power Supply

NEW

MP-310 / MP-320 / MP-510



Cat. No.	MP-100	MINI-300	MINI-500	MP-310	MP-320	MP-510
		 FC CE	 FC CE	 FC CE	 FC CE	 FC CE
Display	N/A	LED		2.4" TFT LCD		
Output voltage / Inc.	50V/100V	10-300V/1V	10-500V/1V	5-300V/1V		5-500V/1V
Output current / Inc.	400mA	10-400mA/1mA	1-400mA/1mA	1-700mA/1mA	10-3000mA/10mA	1-800mA/10mA
Output power / Inc.	40W	60W	120W	150W	300W	
Output type	Constant voltage	Constant voltage/Current		Constant voltage/Current /Power		
Timer	N/A	999 (min) with alarm/Continuous		Constant mode:9999 (min) with alarm/Continuous Programmable mode: 999 (min) with alarm/Continuous		
Rated voltage	100V-120V~ Only			100V-240V~		
Program	N/A	2-steps (V, mA, Ⓢ)		Pre setting; Up to 6-steps (V, mA, W, Ⓢ), 30 programmed files		

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No.19, Ln. 207, Huakang St.
Bade Dist., Taoyuan City 33464
Taiwan
T/ +886-3-3623319
F/ +886-3-3623133

Shanghai Office

Room 612, International business exhibition center,
9300 Hunan Road, Pudong, Shanghai, China
National toll-free No.:400-823-9177
T/ +86-21-50795277
F/ +86-21-50795277

India Office

D. No. 12-13-99, Satguru Apts. Extn.
Street. No. 3, Lane No.1, Tarnaka
Secunderabad – 500 017, India
T/ +91-40-27001515
T/ +91-40-27001586



www.majorsci.com
info@majorsci.com