




3D Rotating Mixer with Adjustable Speed MS-3D-RM

Specification

Model	MS-3D-RM	
Voltage	100~240V, 50/60Hz	
Power	4W	
Setting time	1-999min/0-999s	
Speed	10-40rpm	
Operation angle	20-99°	
Working Mode	Continuous/Timed	
Movement Mode	Rotation, Shaking, Angle	
Application tube (single piece)	50mL	4 wells
	7/10/15mL	5 wells
	5mL	8 wells
	1.5/2mL	24 wells
	0.5mL	11 wells
Machine size	Approx. 240 x 140 x 210mm	
Weight	Approx. 850g	
Certification	CE	



For 0.5 / 50ml tube

For 5 / 7 / 10 / 15ml tube

For 1.5 / 2ml tube

3 different type of tube frames (tube clamps)

Ordering Info

Cat. No.	Product Description
MS-3D-RM-UK	3D Rotating Mixer with Adjustable Speed, UK
MS-3D-RM-EU	3D Rotating Mixer with Adjustable Speed, EU
MS-3D-RM-AU	3D Rotating Mixer with Adjustable Speed, AU
MS-3D-RM-US	3D Rotating Mixer with Adjustable Speed, US

* All images are for reference only, actual products might differ from the pictures above.

* Technical specifications subject to change without notice.

* The tube frame can be different type on the two sides of the cross bar only on the both have same loading. Make sure the samples are in as balanced a position as possible. Please be aware that unbalanced force may cause damage to the machine.

Opional Accessory

Cat. No.	Product Description
MS-3D-RM-001	O-RING, OD14xID7.8mm, MS-3D-RM

Innovative Life Sciences Tools

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Room 612, International business exhibition center,
9300 Hunan Road, Pudong, Shanghai, China
National toll-free No.:400-823-9177
T/ +86-21-50795277
F/ +86-21-50795277

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

India Office

D.No.7-143, 2nd Floor,
St.No.2,Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
T/ +91-40-27001586



*Could be compatibled with MO-MINI.