


3D Mini Tray Shaker

MS-3D-TS series

NEW




MS-3D-TS-001




major science

www.majorsci.com
info@majorsci.com



MS-3D-TS



MS-3D-TS-10


Two rotating angle models are available for different application

Features

- Could be compatible with the MO-MINI.
- Three-dimensional movement is available.
- Two rotating angle models are available for different application

Specification

Name	3D Mini Tray Shaker	
Cat. No.	MS-3D-TS	MS-3D-TS-10
Angle	20°	10°
Movement	Three-Dimensional Movement	
Speed range	0~30rpm	
Speed display	scale (min/max)	
Platform size(WxL)	Approx. 9.21"x6.61"(168x234mm)	
Overall size(WxLxH)	Approx. 6.93"x8.66"x6.5"(176×220×165mm)	
Max Loading	2 kg	
Weight	1.2 kg	
Temperature Range	-10°C~60°C	
Power	20W	
Input voltage	AC100~240V 50/60Hz	
Certification	CE	



*Could be compatible with MO-MINI.

Accessories

Cat. No.	Description
MS-3D-TS-001	optional tray with spring and pad for 3D Tray Shaker, MS-3D-TS.
MS-3D-TS-002	3D Tray Shaker Standard pallet anti-slip mat, rubber, black
MS-3D-TS-003	3D Tray Shaker Optional tray anti-slip mat, rubber, black
MS-3D-TS-004	3D Tray Shaker Optional tray anti-slip mat with spring pull cord, black, 2mm diameter, 350mm circumference.
MS-3D-TS-005	3D Tray Shaker Optional tray anti-slip mat with spring pull cord, black, 2mm diameter, 400mm circumference.

Ordering Info

Cat. No.	Description
MS-3D-TS-CN	3D Tray Shaker with Adjustable Speed, CN power cord
MS-3D-TS-US	3D Tray Shaker with Adjustable Speed, US power cord
MS-3D-TS-UK	3D Tray Shaker with Adjustable Speed, UK power cord
MS-3D-TS-AU	3D Tray Shaker with Adjustable Speed, AU power cord
MS-3D-TS-EU	3D Tray Shaker with Adjustable Speed, EU power cord
MS-3D-TS-IN	3D Tray Shaker with Adjustable Speed, IN power cord
MS-3D-TS	3D Tray Shaker with Adjustable Speed

* All images are for reference only, actual products might differ from the pictures above.
* Technical specifications subject to change without notice