


SafeBlue Electrophoresis System / SafeBlue Imager System



FC CE



Cat. No.	MBE-150	MBE-150-PLUS
Light control	On / Off switch	
Output voltage	10-150V / 1V	
Output current	10-300mA / 1mA	
Output power	30W	
Output type	Constant voltage or constant current	
Timer	999 (min) with alarm / Continuous	



DNA Stains



Midori Green Direct



MIDORI Green Xtra



MIDORI Green Easy



Major Red Plus stain MRP-500



www.majorsci.com
info@majorsci.com

- Innovative - a complete 4-in-1 system (run, view, capture and cut gel) with simple and easy set up on the image capture/analysis software.
- Convenient - run as you see, real time visualization on your DNA/RNA band migration
- Compact - footprint smaller than most UV transilluminator, easy to move and ideal for limited bench top space
- Economical - an inexpensive imager system with digital camera. Multiple accessories are also available to
- increase experimental efficiency

Dual LED Blue/White Light Transilluminator, MBE-200BW

CE

Cat. No.	MBE-200BW
Dimension(mm)(WxLxH)	Approx. 200x200x15.6
Viewing Area(mm)(WxL)	153x153
White Light Wavelength	whole-wavelength
Illuminator Base Design	Flat Bed
Blue Light Wavelength	470nm
Automatic Shutdown	Approx. 6 min



- Combines blue & white full-spectrum LEDs, integrated on opposite sides of a flat-bed design
- Portable size & lightweight
- Real-time observation
- Efficient & early mistake detection
- Ultra-high light uniformity & intensity
- Aluminum alloy casing design
- Low heat dispersion
- Energy saving product (Auto shut-off & LED light)
- Safe & harmless 470 nm blue light wavelength for direct human contact





BluView Transilluminator MBE-200A

CE

- Thinner and more lightweight body
- Easy to carry
- Aluminum alloy casing design
- Low heat dispersion
- Energy saving product
- 470nm harmless blue light for direct human contact

Cat. No.	MBE-200A
Dimension(mm)(WxLxH)	Approx. 200x200x13.9
Viewing Area(mm)(WxL)	153x153
Blue Light Source	15W
Illuminator Base Design	Flat Bed
Blue Light Wavelength	470nm
Automatic Shutdown	Approx. 6 min




BluView Transilluminator MBE-300

CE

- Portable 470nm blue light
- Slim design fits in most of electrophoresis tanks
- Real time observation
- Efficient and early stage mistake detection
- Ultra high light uniformity
- Fit in a majority of different brands of mini and midi size

Cat. No.	MBE-300
Dimension(mm)(WxLxH)	Approx. 86x170x25
Viewing Area(mm)(WxL)	Approx. 112x74.6
Blue Light Source	20W
Illuminator Base Design	Flat Bed
Blue Light Wavelength	470nm
Automatic Shutdown	Approx. 6 min

MBE-300



* All images are for reference only, actual products might differ from the pictures above.
* Technical specifications subject to change without notice

US Office

19959 Sea Gull Way
Saratoga, CA 95070
U.S.A.
T/ +1-408-366-9866
F/ +1-408-446-1107

Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,
Taoyuan City 330041, Taiwan
T/+886-3-3762878
F/+886-3-3761310

Shanghai Office

Building #C1, 880, JiangYang Nan Road,
Baoshan, Shanghai, china 200439
Tel: 13585666474

India Office

D.No.7-143, 2nd Floor,
St.No.2,Nagendra Nagar,
Habsiguda, Hyderabad-500007.
India
T/ +91-40-27001515
F/ +91-40-27001586

