



## Microplate centrifuge

### ▶ MS-MC

- Simple and user-friendly operation.
- Top lid opens to 95° for convenient microplate loading and unloading.
- Compatible with various types of PCR plates (with or without skirt).
- Two operation modes: momentary (inching) and timed run.
- Low noise level during operation.
- Supports rapid, instantaneous spin-down.

Cat. No.	MS-MC
Voltage	AC100~240V, 50/60Hz
Display	Digital
Rotating speed	500~3000 rpm
Timing range	99min 99sec
Max(RCF)	500 xg
Brake time	20s ± 5s
Short run	Yes
Capacity	2× 96-well microplates
Power	45W
Dimension(WxLxH)	220x245x200(mm)
Weight	2.5 kg

\* All images are for reference only, actual products might differ from the pictures above.

\* Technical specifications subject to change without notice.

#### US Office

19959 Sea Gull Way  
Saratoga, CA 95070  
U.S.A.  
T/ +1-408-366-9866  
F/ +1-408-446-1107

#### Taiwan Office

No. 156, Sec. 1, Guoji Rd., Taoyuan Dist.,  
Taoyuan City 330041, Taiwan  
T/+886-3-3762878  
F/+886-3-3761310

#### Shanghai Office

Room 2519, No. 306 Taopu Road,  
Putuo District, Shanghai  
Postal code: 200333  
Cell:+86-13621808204  
T/ +86-21-50795277  
F/ +86-21-50795277

#### India Office

D.No.7-143, 2nd Floor,  
St.No.2,Nagendra Nagar,  
Habsiguda, Hyderabad-500007.  
India  
T/ +91-40-27001515  
T/ +91-40-27001586

